

# Ameda Penguin

## Deluxe Nutritional Warmers

**EN** Instructions for Use 1-15

**ES** Instrucciones de uso 17-32

**FR** Mode d'emploi 33-48



 **Ameda Penguin®**  
Nutritional Warmer

## WELCOME TO THE AMEDA FAMILY

We are excited to help you meet your nutritional warming needs. With the Ameda Deluxe Single-Well or Four-Well Penguin® Nutritional Warmers, you can safely warm feedings and protect vital nutrients.

The Ameda Deluxe Nutritional Warmers work with the patented design Ameda Therma-Liner® which isolates the feeding container from the water while providing an encompassing kangaroo-like warming environment using water to safely transfer heat.

For any questions regarding your Ameda Deluxe Penguin Nutritional Warmer, Call +32 2 304 52 79 or email us at [internationalsales@ameda.com](mailto:internationalsales@ameda.com). Have your serial number for warmers and lot number for Therma-Liners when calling. Please check [www.ameda.com](http://www.ameda.com) for up-to-date hours of operation.

# TABLE OF CONTENTS

<b>1 Warnings And Cautions</b> .....	4-5
<b>2 Intended Use</b> .....	6
<b>3 Product Description</b> .....	6
3.1 Button Description .....	6
<b>4 Set-Up Instructions</b> .....	7
4.1 Preparing The Warmer .....	7
4.2 Customizing Room/Bed/Pod Information.....	7
4.3 Modes .....	7-9
4.4 Therma-Liner Preparation.....	9
<b>5 Inspecting Before Operation</b> .....	10
<b>6 Operating Instructions</b> .....	10
6.1 Warm Refrigerated Or Frozen Feedings .....	10
6.2 Thaw Frozen Milk To Cool Liquid State .....	11
<b>7 Cleaning</b> .....	11
<b>8 Maintenance</b> .....	12
<b>9 Servicing And Spare Parts</b> .....	12
<b>10 Warranty</b> .....	12
<b>11 Technical Data, Storage And Operating</b> .....	13
<b>12 Conditions</b> .....	13
<b>13 Ameda Accessories</b> .....	13



Warning



Temperature Limitation For Transport/Storage



Serial Number



Caution



Humidity Limitation For Transport/Storage



Catalogue Number



Important



Temperature Limitation For Operation



Unit of Measure



Hot Surface



Humidity Limitation For Operation



UL Listed



Grounded



Pressure Limitation For Operation/Transport/Storage



Recyclable  
(Therma-Liners)



Keep Dry



Waste Electrical & Electronic Equipment - the internal components of this product may contain hazardous materials. Dispose according to local or regional waste administration system and regulations. Do not dispose as unsorted municipal waste.



CE marking is a certification mark that indicates conformity with health, safety, and environmental protection standards for products sold within the European Economic Area (EEA)



Manufacturer

# 1. WARNINGS AND CAUTIONS

## IMPORTANT SAFEGUARDS

When using electrical products, the following basic safety precautions should always be followed.

### READ ALL INSTRUCTIONS BEFORE USING THIS PRODUCT

**DANGER** To reduce the risk of electrocution:

- DO NOT place or store where warmer can fall or be pulled into a sink.
- DO NOT place or drop warmer into water or other liquid.
- DO NOT reach for electrical product that has fallen into water. Immediately unplug from electrical outlet.

The **WARNING** symbol identifies important safety instructions that if not followed can lead to injury or damage to the warmer. The warning symbols when used with the following words stand for:



**WARNING** Can lead to serious injury or death.



**CAUTION** Can lead to minor injury.



**IMPORTANT** Can lead to material damage.

## WARNING To avoid fire, electrocution, serious burns or injury.

- **ENSURE** the power cord is easily accessible and removable, in the event of an emergency requiring immediate and total power disconnection.
- **DO NOT** use the equipment in a manner not specified.
- **DO NOT** use AC power source that applies more than 240VAC (nominal) between the supply conductors or between supply conductor and ground. A protective ground connection by way of the grounding conductor in the power cord is required for safe operation.
- **DO NOT** use parts or accessories other than those that are recommended by the manufacturer.
- **USE ONLY** the power cord appropriate for the voltage and plug configuration in your country.
- If the equipment appears damaged or operates abnormally, protection may be impaired. **DO NOT** attempt to continue operation. Have the equipment serviced.

## WARNING To avoid fire, electrocution, serious burns or injury.

- **ALWAYS** operate warmer on a flat, horizontal surface.
- Close supervision is necessary if product is used near children or infants.
- Use product for its intended use **ONLY**, as described in this manual.
- Prior to use, always inspect power cord for damage or exposed wires. If damage and/or exposed wires are found, **DO NOT** use and call us at +32 2 304 52 79.
- **NEVER** operate electrical warmer
  - if it has a damaged cord or plug
  - if it is not working properly
  - if it has been dropped, damaged, or becomes wet without being properly and fully assembled
- **ALWAYS** keep cord away from heated surface.
- **ALWAYS** ensure the power cord is disconnected before cleaning or servicing.
- **ALWAYS** ensure the voltage of the product is compatible to the power source.
- **ALWAYS** ensure the power cord is inserted in a socket-outlet with a protective earth contact for grounding. Intentional AC interruption is prohibited.
- **ALWAYS** ensure servicing is performed by qualified service personnel only using the Manufacturer's Service Manual. There are no user-serviceable parts.

## CAUTION Can lead to minor injury.

- **DO NOT** use Therma-Liner
  - for more than one shift, not to exceed 12 hours
  - on more than one patient, baby, or child

## IMPORTANT Can lead to material damage.

- **NEVER** put electrical product in water or a sterilizer, as you can cause permanent damage to the nutritional warmer.
- **DO NOT** use abrasive cleaners/detergents when cleaning the nutritional warming unit or accessories. CaviWipe™ or SaniWipe™ are recommended for surface cleaning.
- **DO NOT** use sharp objects on the device buttons or wells.
- **DO NOT** use Therma-Liner for more than one shift or 12 hours.
- **DO NOT** overfill or under-fill Therma-Liner.
- If servicing is required, ensure servicing is performed by qualified service personnel.



**SAVE THESE INSTRUCTIONS**

## 2. INTENDED USE










The Ameda Penguin® Nutritional Warmer is intended to warm feedings to a range that includes body temperature.

## 3. PRODUCT DESCRIPTION Penguin Nutritional Warmer


Your Penguin Nutritional Warmer is a feeding warming device and is customized with settings and warming profiles according to your unique facility needs. The nutritional warmer operates with gentle vibration which helps preserve and suspend feeding nutrients and promote even warming.

### 3.1 BUTTON DESCRIPTION

#### Navigation Buttons

Top		Scroll Up		Scroll Up	
Middle		Select		Go Back	
Bottom		Scroll Down		Scroll Down	

#### Beak Buttons




## 4. SET-UP INSTRUCTIONS

### 4.1 PREPARING THE WARMER

1. **PLACE** warmer onto a stable, flat, horizontal surface.
2. **PLUG** the power cord into the power source.
3. **CONNECT** the power cord into the back of the warmer device.

### 4.2 CUSTOMIZING ROOM/BED/POD INFORMATION


1. **FLIP** the rear power switch into the “**ON**” position.
2. **PRESS** and **HOLD** the bottom navigation button.
3. **FOLLOW** the on-screen instructions to ENTER up to 5 spaces for Room, up to 5 spaces for Bed and up to 5 spaces for Pod:
  - +abc/123
  - =return to previous entry
  - -abc/123
  - Enter button=advance to next entry
  - On/Off button=edits complete, return to Standby mode




### 4.3 MODES

#### STANDBY MODE

The Penguin Nutritional Warmer will be in Standby mode and cycle through the following screens whenever powered on and not warming:




## WARMING MODES



### Fridge to Feed

**Feed Temp** range of 90-100°F/32-38°C


**Average Time** 5-12 minutes depending on volume



### Frozen to Feed

**Feed Temp** range of 90-100°F/32-38°C

**Average Time** 10-18 minutes depending on volume



### Thaw

**After Thaw Temp** range of 38-45°F/3-7°C

**Time** 15-25 minutes depending on volume

**READY TO FEED MODE - HOLD TIME** (Fridge to Feed and Frozen to Feed only) The feeding is ready when the screen background turns green with the following words displayed:




### REMAINING HOLD TIME

Displays a countdown of the remaining hold time. The warmer emits a ready tone at the beginning of the hold time and periodically throughout. The warmer will keep the feeding in a safe temperature range for 30 minutes (until turned off, returned to Standby mode, or extends to Overtime mode).



## THAW COMPLETE MODE



When the Thaw feature is used, the feeding should be removed immediately upon completion to prevent over-warming. The warmer emits a ready tone. The heater immediately shuts off, and milk must be refrigerated. The screen background will turn green and the phrase shown will appear on the screen.

**WARNING** Milk must be refrigerated immediately.

## OVERTIME MODE Fridge to Feed and Frozen to Feed only



After 30 minutes in Ready to Feed HOLD TIME, the warmer will enter Overtime. Warming and vibration will turn off. The screen background will turn red and the word "Overtime" will appear in the middle of the screen.

The Overtime screen will remain on until the unit is turned off or is returned to Standby by PUSHING the On/Off button or Enter button.

**Note: As the manufacturer, we recommend not using the feeding after the 30 minute hold time, but that feeding should be handled according to the hospital policy.**

### DISPLAYED ERRORS

If the Penguin Nutritional Warmer encounters an error, the error will be displayed on a red screen with a warning tone. If this should occur, discontinue use of the warmer and call Ameda CustomerCare at 1.877.99.AMEDA (26332).

## 4.4 THERMA-LINER PREPARATION

1. At the beginning of a shift, take out one Therma-Liner per user.
2. Using a permanent marker, enter the following information onto the Therma-Liner:
  - a. Patient Information
  - b. Date/Time
  - c. Information
3. Open the water compartment (blue zipper) and fill to water line using room temperature water.
4. Pinch one end of the blue zipper between your thumb and forefinger and slide down the length of the zipper to close the water compartment.
5. Repeat to ensure a proper seal. This compartment will remain sealed for the duration of the 12 hour shift.

**Note: Use care when removing and re-inserting milk container to prevent damage to interior of water compartment.**

## 5. INSPECTING BEFORE OPERATION

Inspect the product for damage. Do not use if damaged.  
Call Ameda CustomerCare at 1.877.99.AMEDA (26332).

## 6. OPERATING INSTRUCTIONS

1. Flip the rear power switch into the 'ON' position.
2. Press On/Off Button to begin.
3. Scroll to the warming mode of interest:
  - a. Room Temperature/Fridge to Feed
  - b. Frozen to Feed
  - c. Thaw
4. Advance by pressing Enter or Select (green check mark)\*.

*\*Note that if a wrong entry was made, you may press the Go Back selection button to go back a screen.*

### 6.1 WARM REFRIGERATED OR FROZEN FEEDINGS

1. Scroll to choose the container type and size of interest.
2. Press Enter.
3. Only when ready to begin warming the feeding, place the feeding container into the inner compartment (orange zipper) of the Therma-Liner and gently push it to the bottom. Press firmly along the entire orange zipper to encompass the container as much as possible. Use care when removing and re-inserting feeding containers to prevent damage to the interior water compartment.
4. Place the prepared Therma-Liner in the well with the sealed opening at the top.
5. Press Enter\*.
6. The device will begin warming, vibrating, and alternating screens between the profile selected and elapsed time.
7. When the feeding has been warmed to a range that includes body temperature, the device will enter Ready to Feed mode.
8. Please remove the feeding from the Therma-Liner and set the Therma-Liner aside until next use during the shift.
9. Press On/Off to return the Penguin to Standby mode.

*\*Note that if a wrong entry was made, you may press the Cancel selection button to go back a screen.*

## 6.2 THAW FROZEN MILK TO COOL LIQUID STATE

*Caution: the Penguin Nutritional Warmer is not recommended for thawing less than 15 mL of breast milk in storage bags!*

1. Scroll to choose the container type and size of interest.
2. Press Enter.
3. Only when ready to begin warming the feeding, place the feeding container into the inner compartment (orange zipper) of the Therma-Liner and gently push it to the bottom. Press firmly along the entire orange zipper to encompass the container as much as possible. Use care when removing and re-inserting feeding containers to prevent damage to the interior water compartment.
4. Place the prepared Therma-Liner in the well with the sealed opening at the top.
5. Press Enter\*.
6. The device will begin warming, vibrating, and alternating screens between the profile selected and elapsed time.
7. When the feeding has been warmed to a liquid cool state, the device will enter Thaw Complete mode.
8. Please immediately remove the feeding\*\* from the Therma-Liner and set the Therma-Liner aside until next use during the shift.
9. Press On/Off to return the warmer to Standby mode.

*\*Note that if a wrong entry was made, you may press the Cancel selection button to go back a screen.*

*\*\*If there are still some frozen particles in the feeding, let it sit until the remainder thaws, or rotate the container gently to thaw the rest of the way. It is better to have a few remaining ice crystals than to overshoot the refrigerated range. Any portion that is meant for immediate feeding may be placed in a bottle or syringe for further warming. Feeding for the future should be placed immediately in the refrigerator for up to 24 hours. After 24 hours, it should be discarded.*

## 7. CLEANING

1. Disconnect Power Cord from the wall receptacle.
2. Prepare reagent wipes (e.g. CaviWipe™ or SaniWipe™) containing Isopropanol, Ethylene Glycol Monobutyl Ether (2-Butoxyethanol), Diisobutylphenoxyethoxyethylidimethyl benzylammonium chloride.
3. Using the prepared reagent wipe, clean and remove all surface debris from the device, including the power cord.
4. Repeat as deemed necessary using a newly prepared reagent wipe each time.
5. After the last cleaning wipe is used, allow device to air dry for a minimum of 10 minutes or dry with a soft cloth before returning the device to service.

## 8. MAINTENANCE

The Penguin Nutritional Warmer requires no preventative maintenance.

## 9. SERVICING AND SPARE PARTS

All warmer service and repair can be accomplished by contacting Ameda. CustomerCare will issue a Return Material Authorization (RMA) or provide your facility's biomedical department with parts needed to repair the warmer and return it to operation.

## 10. WARRANTY

### **THE AMEDA WARRANTY IS NONTRANSFERABLE. THE AMEDA PENGUIN NUTRITIONAL WARMER IS WARRANTED TO THE ORIGINAL OWNER ONLY.**

From the date of the purchase of this product, Ameda, Inc. ("Ameda") warrants the Product against defects in material or workmanship for a period of two years. Ameda's sole obligation under this limited warranty shall be to repair or replace, at Ameda's option, any Product that is determined to be defective by Ameda and determined to be covered by this express limited warranty. Repair or replacement under this limited warranty is the sole and exclusive remedy of the First User. PROOF OF PURCHASE IN THE FORM OF A RECEIPTED INVOICE OR BILL OF SALE EVIDENCING THAT THE PRODUCT IS WITHIN THE WARRANTY PERIOD MUST BE PRESENTED TO OBTAIN WARRANTY SERVICE. This limited warranty is extended by Ameda ONLY to the First User and is not assignable or transferable. For warranty service contact Ameda at 1.877.99.AMEDA (26332).

EXCEPT TO THE EXTENT PROHIBITED BY APPLICABLE LAW, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE ON THE PRODUCT IS HEREBY DISCLAIMED. AMEDA SHALL NOT BE LIABLE FOR ANY INCIDENTAL, CONSEQUENTIAL, SPECIAL OR PUNITIVE DAMAGES OF ANY KIND FOR BREACH OF THE EXPRESS LIMITED WARRANTY ON THE PRODUCT OR ANY WARRANTY IMPLIED BY OPERATION OF LAW. OTHER THAN THE EXPRESS LIMITED WARRANTY SET FORTH ABOVE, THERE ARE NO OTHER WARRANTIES THAT ACCOMPANY THIS PRODUCT AND ANY ORAL, WRITTEN OR OTHER REPRESENTATION OF ANY KIND IS HEREBY DISCLAIMED.

The user/owner of this equipment shall have sole responsibility and liability for any injury to person or damage to property resulting from: 1) Operation not in accordance with supplied operating instructions, 2) Maintenance not in accordance with authorized maintenance service instructions, 3) Service by anyone other than manufacturer authorized qualified service personnel, 4) Modification of the equipment or accessories and 5) Use of damaged or unauthorized components and accessories.

## 11. TECHNICAL DATA, STORAGE AND OPERATING CONDITIONS

### **Penguin Single Well Nutritional Warmer**

Power Rating: 220-240VAC, 50/60 Hz, 1.25A

Relative Humidity: 80% up to 31°C to 50% to 40°C / 80% up to 87.8°F to 50% up to 104°F

Altitude: 2000 meters

Operating Temperature: 5 to 40°C/41 to 104°F

Storage Temperature: -29 to 70°C/-20 to 158°F

Line Fuse: 2 each, type "T" 1.25A, 250VAC

Mass: 3.96 kg

Size: 21.8 cm H X 19.4 cm W X 21.9 cm D

Pollution degree 2

Installation category H

Indoor use only

## 12. AMEDA ACCESSORIES

You may also order additional warmer accessories: Therma-Liner replacements or Perche Shelf through your sales specialist.



Therma-Liner



Perche Shelf

# 13. EMC INFORMATION

THIS **ELECTRICAL EQUIPMENT** needs special precautions regarding EMC and needs to be installed and put into service according to the EMC provided in these Instructions for use. Portable and mobile RF communications equipment can affect this **ELECTRICAL EQUIPMENT**. **Warning:** Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Penguin Nutritional Warmer, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result. **Warning:** Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation. **Warning:** Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.


Guidance and manufacturer's declaration – electromagnetic emissions		
The Ameda Penguin Warmer is intended for use in the electromagnetic environment specified below.		The customer or the user of the Ameda Penguin Warmer should assure that it is used in such an environment.
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The Ameda Penguin Warmer uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration – electromagnetic immunity			
The Ameda Penguin Warmer is intended for use in the electromagnetic environment specified below.		The customer or the user of the Ameda Penguin Warmer should assure that it is used in such an environment.	
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/ burst IEC 61000-4-4	± 2 kV for AC line; 100 kHz repetition frequency	± 2 kV for AC line; 100 kHz repetition frequency	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 0,5 kV, ± 1kV line-to-line; ± 0,5 kV, ± 1kV, ± 2kV line-to-ground	± 0,5 kV, ± 1kV line-to-line; ± 0,5 kV, ± 1kV, ± 2kV line-to-ground	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<ul style="list-style-type: none"> <li>● 0% Ut; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°</li> <li>● 0% Ut; 1 cycle and 70% Ut; 25/30 cycles Single phase: at 0°</li> <li>● 0% Ut; 250/300 cycles</li> </ul>	<ul style="list-style-type: none"> <li>● 0% Ut; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°</li> <li>● 0% Ut; 1 cycle and 70% Ut; 25 cycles At 0°</li> <li>● 0% Ut; 250 cycles</li> </ul>	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Ameda Penguin Warmer requires continued operation during power mains interruptions, it is recommended that the Ameda Penguin Warmer be powered from an uninterruptible power supply.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m 50 Hz or 60 Hz	30 A/m 50 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
<b>NOTE</b> Ut is the a.c. mains voltage prior to application of the test level.			

## Guidance and manufacturer's declaration – electromagnetic immunity

The Ameda Penguin Warmer is intended for use in the electromagnetic environment specified below.

The customer or the user of the Ameda Penguin Warmer should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 V 0,15 MHz – 80 MHz 80% AM at 1 kHz	3 V 0,15 MHz – 80 MHz 80% AM at 1 kHz	Portable and mobile RF communications equipment should be used no closer to any part of the Ameda Penguin Warmer, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1,2\sqrt{P}$
	6V in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80% AM at 1 kHz	6V in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80% AM at 1 kHz	
Radiated RF IEC 61000-4-3	10 V/m 80 MHz – 2,7 GHz 80% AM at 1 kHz	10 V/m 80 MHz – 2,7 GHz 80% AM at 1 kHz	$d = 1,2\sqrt{P}$ 80 MHz to 800 MHz $d = 2,3\sqrt{P}$ 800 MHz to 2,7 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, <sup>a</sup> should be less than the compliance level in each frequency range. <sup>b</sup>  Interference may occur in the vicinity of equipment marked with the following symbol: 
	9 – 28 V/m At IEC 60601-1-2: 2014 Table 9 frequencies	9 – 28 V/m At IEC 60601-1-2: 2014 Table 9 frequencies	

**NOTE 1** At 80 MHz and 800 MHz, the higher frequency range applies.

**NOTE 2** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Field strength from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the

electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Ameda Penguin Warmer is used exceeds the applicable RF compliance level above, the Ameda Penguin Warmer should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Ameda Penguin Warmer.

b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

## Recommended separation distances between portable and mobile RF communications equipment and the Ameda Penguin Warmer

The Ameda Penguin Warmer is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Ameda Penguin Warmer can help prevent electromagnetic interference by maintaining a minimum

distance between portable and mobile RF communications equipment (transmitters) and the Ameda Penguin Warmer as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = 1,2\sqrt{P}$	80 MHz to 800 MHz $d = 1,2\sqrt{P}$	800 MHz to 2.7 GHz $d = 2,3\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

**NOTE 1** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

**NOTE 2** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.



**Distributed in USA by:**

Ameda, Inc.  
485 Half Day Road  
Buffalo Grove, IL 60089  
1.866.99.AMEDA (26332)  
www.ameda.com

**Distributed in Australia by:**


MidMed Pty Ltd.  
Unit 4/62 Borthwick Avenue  
Murarrie QLD 417261 (7)  
3348 9155  
www.midmed.com.au

Distribution in other countries: For product information or feedback, call your local distributor or location where you purchased the product. For a listing of distributors in your country, please visit [www.ameda.com](http://www.ameda.com).


Distribution dans d'autres pays : Pour plus de renseignements ou d'avis sur le produit, communiquez avec votre distributeur local ou l'endroit où vous l'avez acheté. Pour obtenir une liste de distributeurs dans votre pays, veuillez visiter [www.ameda.com](http://www.ameda.com).


Distribución en otros países: para obtener información o comentarios sobre el producto, llame a su distribuidor local o al lugar donde compró el producto. Para obtener una lista de los distribuidores en su país, visite [www.ameda.com](http://www.ameda.com).


**Transport/Storage**  
**Transporte/almacenamiento**  
**Transport/Rangement**



**Operation**  
**Funcionamiento**  
**Fonctionnement**



 **Manufactured for:**  
**Ameda, Inc.**  
485 Half Day Road,  
Buffalo Grove, IL 60089 USA



1500165 0819