



HygieniKit

MILK COLLECTION SYSTEMS



CustomFit Options | Breast pump “flange fit” is vital to a mother’s comfort and milk flow. But achieving a good fit with other brands costs time and money. Now there is no need to carry several kit configurations to achieve optimal fit. Ameda HygieniKits are available with multiple flange sizes to make finding a good fit easier—at first and later as a mother’s body changes.

Easy Care | Why use a breast pump kit that requires tubing inspection after every use? The Ameda HygieniKit Milk Collection system features a patented diaphragm that provides a solid barrier between the milk and the tubing. Moisture and milk cannot enter tubing, simplifying pump use and care and reducing risk of infection in compromised babies.

Compatibility | With some brands, when a mother changes pumps, this means new kit pieces or configurations are needed. Not with Ameda. All Ameda HygieniKits attach the same way to all Ameda electric breast pumps.

Key Features

- Easy Care & Cleaning
- Easy to Use—Sterile, Ready to Use
- Less Inventory to Meet All Pumping Needs
- Reduced Risk to Compromised Babies
- CustomFit Flange Options
- Converts Easily to Manual Pump
- Pump Into Milk Storage Bag or Bottles
- Single, Dual and Manual Kits Available

Adaptability & Simplicity



MOM INSPIRED. HOSPITAL TRUSTED.

HYGIENIKITS with CustomFit Flanges

(Dual and Single HygieniKits with 30.5 mm, 28.5 mm and 25.0 mm Flange Sizes)



17449 Dual HygieniKit with CustomFit Flanges and One-Hand Manual Breast Pump Adapter



17448 Single HygieniKit with CustomFit Flanges and One-Hand Manual Breast Pump Adapter

All
Ameda HygieniKits
are without BPA

CUSTOM BREAST FLANGES with Inserts

(36.0 mm, 32.5 mm, 30.5 mm and 28.5 mm Flange Sizes)



17358 XX-Large Custom Breast Flanges with X-Large Inserts (36.0 mm/32.5 mm) (Non-sterile)



17148P Large Custom Breast Flanges with Medium Insert (30.5 mm/28.5 mm) (Non-sterile)

Onekit
Manyoptions

Every Ameda HygieniKit is designed to fit All Ameda Breast Pumps

DUAL HYGIENIKIT MILK COLLECTION SYSTEMS

(25.0 mm Standard Breast Flange)



17152 Dual HygieniKit Milk Collection System with One-Hand Manual Breast Pump Adapter



17155 Dual HygieniKit Milk Collection System



17153 Dual HygieniKit Milk Collection System with One-Hand Manual Breast Pump Adapter and Flexishield Areola Stimulators

SINGLE HYGIENIKIT MILK COLLECTION SYSTEMS

(25.0 mm Standard Breast Flange)



17150 Single HygieniKit Milk Collection System with One-Hand Manual Breast Pump Adapter



17151 Single HygieniKit Milk Collection System

CustomFit Flange System™

Every pumping mother needs a good flange fit for greater comfort and better milk flow.

- Flange fit is based on the size of the flange opening —or nipple tunnel— the breast is drawn into during pumping.
- Flange fit can change with birth, breastfeeding and pumping.

To help every mother achieve her best fit, Ameda's CustomFit Flange System™ offers seven flange sizes – ranging from 21.0 mm to 36.0 mm – that work with the Ameda HygieniKits.

For your convenience, some Ameda HygieniKits include three flange sizes to make finding a good fit even easier (see top left).



For more information about the Ameda HygieniKit Milk Collection System call 1.877.99.AMEDA (26332) or visit ameda.com

Each HygieniKit shown may not be available in all countries. Check with your local distributor for details.

Distributed in USA by:
Evenflo Company, Inc.
1801 Commerce Drive
Piqua, Ohio 45356
1.877.99.AMEDA (26332)
www.ameda.com

Distributed in Canada by:
Mothers Choice Products Ltd
2133-11871 Horseshoe Way
Richmond, BC V7A 5H5
1.800.604.6225
www.motherschoiceproducts.com

Distributed in Australia by:
Advantage Health Care
34 Percy Street
Mordialloc VIC 3195
61 (3) 9580 9288
www.ameda.com.au

Ameda, CustomFit Flange System, HygieniKit and associated logos are trademarks of Evenflo Company, Inc.
©2008 Evenflo Company, Inc.
26601006-1108